

## Product Name

Dental Polishing Burs

## Structure and Composition

This rotary instrument consists of a metal shank and a working end. The metal shank connects to the dental handpiece, and the working end is shaped as a cylinder, cup, or disc. It is supplied non-sterile.

## Intended Use

For grinding and polishing restorations, removing excess material, as well as cleaning and polishing tooth surfaces.

## Instructions for Use

Use with a dental handpiece. Secure the product onto the handpiece and apply directly to the surface of dentures or oral restorations at an appropriate angle for grinding and polishing.

## Other Operating Precautions

1. Use in accordance with the handpiece manufacturer's instructions. Ensure the product is fully inserted into the handpiece; do not use if incompletely seated.
2. Start at a speed below 2000 r/min and confirm there is no vibration before use.
3. Operate at the recommended speed specified on the outer packaging. The maximum speed shall not exceed 20,000 r/min. For products with smaller/larger diameters or longer lengths, use at a 45° angle. Do not apply excessive force, as this may cause breakage.
4. Disinfectant sterilization may cause oxidation and rust on the surface. Use steam sterilization at 135°C. The product packaging box must not be used as a sterilization container.
5. For intraoral use: single-use only. Repeated sterilization may reduce the bonding strength of the attachment. The product may be reused by technicians in a laboratory setting.
6. Do not use deformed or damaged products.
7. Wear a mask and protective eyewear during use.
8. Do not use for purposes other than those specified in the intended use or on the outer packaging.
9. This product is intended for use by dental professionals only.

## Storage Conditions

Store away from direct sunlight, moisture, and corrosive environments.

## Transportation Conditions

Avoid squeezing and impact during transportation and storage.

## Date of Manufacture

See product label.

## Service Life

10 years.



**AT&M Biomaterials Co., Ltd.**

Building 1, Zhongguancun Science and Technology Park at Pinggu, Courtyard No. 46,

Pingwang Street, Pinggu District, Beijing, China

Tel: +86-10-69778208

Email: [info@atmbio.com](mailto:info@atmbio.com)

Website: [www.atmbio.com](http://www.atmbio.com)